

1-800-548-2362

International (912) 236-0000



www.innomed.net

(912) 236-7766 Fax

©2000 Innomed, Inc., All Rights Reserved.

**Innomed**  
Instrument Division  
103 Estus Drive, Savannah, GA 31404



FREE TRIAL ON MOST INSTRUMENTS  
CALL FOR MORE INFORMATION

## Kirschenbaum Acromioplasty Retractor

Designed to fit under the posterior edge of the acromion and lever the humeral head down out of the way. Helps to protect both the posterior aspect of the shoulder and the articular surface of the humeral head during open acromioplasty and rotator cuff surgery.

*Designed by*  
I. Kirschenbaum, M.D.

Product No: 5840

