


Arm Board

A radiolucent arm board designed for all upper extremity surgery

- ▶ Securely slides under the table pad
- ▶ Patients body weight holds arm board in place



PRODUCT NO:

2419

Overall Length: 42" (106,7 cm)
Maximum Width: 19" (48,25 cm)
Minimum Width: 10" (25,4 cm)
Thickness: .5" (1,27 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

www.innomed.net



INNOMED 



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362