


Articulated Caliper with Ruler


Designed by Vincent Y. Ng, MD

A highly precise (within 1 mm) device designed for measuring distances between two points – can be used even if there are intervening structures like soft tissue or bone, and in situations where a straight ruler will not work

Examples of use include measuring limb length in total hip arthroplasty, confirming length in megaprosthesis knee replacements, and assessing dimensions of allografts.



Articulated Caliper
#2026-01



Ruler #2026-02

PRODUCT NO.'S:	
2026-00	[Articulated Caliper with Ruler]
Set Includes/ Available Separately:	
2026-01	[Articulated Caliper Only]
Overall Length (unfolded): 15.25" (38,8 cm)	
Dimensions Triangle Folded: 4" x 4.25" (10,2 x 10,8 cm)	
2026-02	[Ruler Only]
Overall Length: 9" (22,9 cm)	
Width: .79" (2 cm)	



Reference Measurement Statement (Non-Measuring Function)

This device is not intended to provide quantitative measurements and does not incorporate a measuring function as defined under Regulation (EU) 2017/745 (MDR). Any scales, markings, or visual indicators present on the device are provided solely for reference and orientation purposes to support general use. These indicators are not calibrated against traceable standards and have not been validated for accuracy, precision, or repeatability. Accordingly, values derived from these markings shall be considered approximate only and must not be used where accurate or traceable measurements are required for clinical decision-making, diagnosis, or other applications requiring defined metrological performance. Users shall rely only on devices specifically intended and validated for measurement when quantitative accuracy is required.

www.innomed.net



INNOMED

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362