

Bacastow Glenoid Retractor

Designed by David Bocastow, MD

Designed for glenoid exposure, particularly for reverse shoulder replacement applications, where it is important to get inferiorly


Allows visualization and direct access to the glenosphere base plate through a deltopectoral incision with intact pectoralis major insertion.

PRODUCT NO'S:

1897-L [Left]
Overall Length: 11.75" (29,8 cm)

1897-R [Right]
Overall Length: 11.75" (29,8 cm)

PROUDLY
MADE IN
THE USA



ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Scan to
Launch Our
Website



Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362