


Blount Retractor with Small Handle

Designed by Ronald Romanelli, MD

A blount retractor with a lightweight ergonomic handle designed for tissue retraction and closure assistance in knee, shoulder, and hip arthroplasty



PRODUCT NO:

4852

Overall Length: 9.375" (25,1 cm)

Handle Length: 4.625" (11,7 cm)

Blade Depth: 1.5" (3,8 cm)

Blade Width at Widest: .325" (8,3 mm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362