





Burkhead Glenoid Retractor

Designed by Wayne "Buzz" Burkhead, Jr, MD

The retractor bar presses against the glenoid while the end of the retractor puts pressure on the posterior capsule

PRODUCT NO'S:
5839 [Large] Overall Length: 9.125" (23,2 cm) Blade Width at End: 1.5" (3,8 cm)
5839-SM [Small] Overall Length: 8.75" (22,2 cm) Blade Width at End: 1" (2,54 cm)

PROUDLY
MADE
IN THE
USA



ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362