

Burkhead Reversible TSA/RSA Retractor

Designed by Wayne "Buzz" Burkhead, Jr, MD

Unique shape, angles and double pronged end serves to push the posterior capsule, and the humerus, away from the glenoid to allow preparation of the glenoid and implantation of component(s) without having to remove the retractor



PRODUCT NO:

5839-01

Overall Length: 9.125" (23,2 cm)

Blade Width at End: 1.5" (3,8 cm)

**PROUDLY
MADE
IN THE
USA**



ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362