


# Calibrated Ortho Spreader with Slotted Tips

Designed by Jason Bariteau, MD

*A lamina spreader with a very thin closed profile, designed to enable distraction in tight joint spaces like the subtalar and talonavicular joints*




**PRODUCT NO:**

**1841**

Overall Length: 6.75" (17,1 cm)  
Prong Length: .5" (12,7 mm)  
Calibrations: 10 mm to 35 mm



*Locations for use of Calibrated  
Ortho Spreader with Slotted Tips*



[www.innomed.net](http://www.innomed.net)



ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

**INNOMED** 



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
info@innomed-europe.com

**1.800.548.2362**

© 2025 Innomed, Inc.