



Chandler Retractor

Used for retracting tissue away from the bone, and helpful for posterior exposure of the tibia in MIS surgery

Allows the surgeon to retract soft tissue away from the bone, and can be used for hip and knee surgery. The handle is contoured away from the field of view and working area. Available in three blade sizes: 5/8", 3/4" and 1".

The OrthoLucent™ version is made of a strong, lightweight carbon fiber PEEK composite material, which is completely radiolucent, helps to prevent from marring component surfaces, and can be steam sterilized.

PRODUCT NO'S:

3220-01 [5/8"] Overall Length: 9.875" (25.1 cm) Blade Width: 5/8" (1.6 cm)	3220-04 [1"] Overall Length: 9.875" (25.1 cm) Blade Width: 1" (2.54 cm)
3220-02 [3/4"] Overall Length: 9.875" (25.1 cm) Blade Width: 3/4" (1.9 cm)	3220-02R* [OrthoLucent™ 3/4"] Overall Length: 9.875" (25.1 cm) Blade Width: 3/4" (1.9 cm)

**PROUDLY
MADE
IN THE
USA**

*
**MADE EXCLUSIVELY
FOR INNOMED IN
SWITZERLAND**



ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362