



Chandran Distal Biceps Tissue Protector


Designed by Rama E. Chandran, MD

Designed to help protect tissue and expose the radial tuberosity during distal biceps tendon repair, and also useful to help expose the humerus during proximal subpectoral biceps repair

Using downward pressure, the teeth help to engage bone to keep the protector in place.


PRODUCT NO:

3224
Overall Length: 8" (20,3 cm)
Handle Length: 5.5" (14 cm)
Blade Depth: 1.7" (4,3 cm)
Blade Width: .65" (1,7 cm)



Helps to protect tissue and expose the radial tuberosity during distal biceps tendon repair.

Also useful to help expose the humerus during proximal subpectoral biceps repair.



www.innomed.net

FREE TRIAL ON MOST INSTRUMENTS



ISO 13485:2016



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.