

1-800-548-2362

International (912) 236-0000



www.innomed.net

(912) 236-7766 Fax

©2006 Innomed, Inc., All Rights Reserved.

**Innomed**

103 Estus Drive, Savannah, GA 31404  
E-Mail: info@innomed.net



FREE TRIAL ON MOST INSTRUMENTS  
CALL FOR MORE INFORMATION


## Chuck Key Handle

The Chuck Key Handle is designed to snap onto a standard chuck key giving better leverage during tightening of a chuck. It also helps keep a chuck key from slipping or being dropped during surgery.

Product No: 5560  
Chuck Key Not Included



*Snaps onto a standard chuck key for better leverage*



## Large Handle Chuck Key

Designed to allow a chuck to be tightened and untightened easily. Available with a 1/4" chuck key on a 7" shaft.

Product No's: 5517-01



*For easy tightening/untightening of a chuck*

