

Cobb Elevators

Ultra hard titanium nitride coating helps to extend blade life by increasing surface hardness, prolonging sharpness, and resisting chemicals and corrosion.



Two Sizes Available
With or Without Teeth

PRODUCT NO'S:	
WITH TEETH	WITHOUT TEETH
3432 [1/2" with Teeth] Overall Length: 11" (27,9 cm) Blade Width: 1/2" (13 mm)	3436 [1/2" without Teeth] Overall Length: 11" (27,9 cm) Blade Width: 1/2" (13 mm)
3434 [1" with Teeth] Overall Length: 11" (27,9 cm) Blade Width: 1" (25,4 mm)	3438 [1" without Teeth] Overall Length: 11" (27,9 cm) Blade Width: 1" (25,4 mm)



www.innomed.net

FREE TRIAL ON MOST INSTRUMENTS

INNOMED

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.