


Connelly Leg Positioning Set


Designed by Camille Connelly, MD

A clamp free leg positioning board for partial or total knee replacement designed to go UNDER the drapes, with a sterilizable footpiece for easy positioning



- ▶ The board can be stabilized on the OR table with Velcro or tape
- ▶ Set includes the footpiece
- ▶ Two (2) Sterile Pads/Wraps are included with each new purchase (See below for Pad/Wrap information)

PRODUCT NO'S:
2626-00 [Connelly Leg Positioning Set]
SET INCLUDES / Available Individually:
2626-01 [Connelly Leg Positioning Board] Overall Length: 30" (76,2 cm) Overall Width: 10" (25,4 cm) Height: 1.25" (3,2 cm)
2626-FP [Connelly Leg Positioner Footpiece] Length: 11.9" (30 cm) Height: 7.9" (20 cm) Width: 6" (15,2 cm)
2629-L [Sterile Pad/Wrap] Two included in set; product number is for one only
Replacement Parts:
2629-00 [Case of 10 Sterile Pads/Wraps]



OPTIONAL PRODUCT NO:
8840 [Leg Stabilizer]



Leg Positioner Sterile Protective Pad & Wrap

Disposable, latex-free sterile foam pad and cohesive wrap helps protect patient from pressure sores, abrasions and possible neurological impairment while securing foot into the boot.

PRODUCT NO'S:
2629-00 [Case of 10 Sets]
2629-L [1 Set]



Each set includes 1 Pad and 1 Wrap

www.innomed.net



FREE TRIAL ON MOST INSTRUMENTS

INNOMED

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2026 Innomed, Inc.