




Dillon Curved Glenoid Retractors

Designed by Mark Dillon, MD

Designed to help bring the retractor away from the glenoid in reverse TSA procedures, helping allow for easier reaming of the glenoid, while the surgeon can also use the curved design at the end to help guide the glenosphere onto the baseplate, making insertion significantly easier



Standard
#4278

Narrow
#4278-N

PRODUCT NO'S:	
4278	[Standard]
Overall Length: 10.875" (27,6 cm)	
Handle Width: 4.33" (11 cm)	
Blade Width: 1.5" (3,8 cm)	
4278-N	[Narrow]
Overall Length: 10.125" (25,7 cm)	
Handle Width: 4.33" (11 cm)	
Blade Width: 1.45" (3,7 cm)	



www.innomed.net



FREE TRIAL ON MOST INSTRUMENTS

INNOMED 

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362