



Sweed Dissecting Scissors

Designed by Tamer Sweed, FRCS (Orth)

Designed with a blunt, flat bar fixed to the lower limb, the scissors also act as a dissector to protect underlying vital structures



PRODUCT NO:

3081

Overall Length: 6.625" (16.8 cm)

Bottom Pad: 16 mm x 6 mm

Pad Extension Beyond Scissor: 6 mm

MADE EXCLUSIVELY
FOR INNOMED IN
GERMANY

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED 



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362