


Ditmars Carpal Tunnel Release Set

Designed by Donald M. Ditmars Jr., MD

Designed to help retract and provide access for carpal tunnel release operations


Large Release Retractor

Retractor for carpal tunnel release operations
7,5 mm inside tube diameter



Small Release Retractor

Small retractor for initial use in carpal tunnel release operations
4 mm inside tube diameter



Straight Carpal Tunnel Probe

Probe to act as dilator to insert small retractor for carpal tunnel release operation



PRODUCT NO'S:

1132-00 [Carpal Tunnel Release Set with Case]

Also Available Individually:

1132-01 [Large Curved Release Retractor]

Overall Length: 5" (12,8 cm)

Handle Length: 3" (7,6 cm)

Inside Tube Diameter: 7,5 mm

1132-02 [Small Curved Release Retractor]

Overall Length: 4.75" (12 cm)

Handle Length: 3" (7,6 cm)

Inside Tube Diameter: 4 mm

1132-03 [Straight Carpal Tunnel Probe]

Overall Length: 7.5" (19,1 cm)

Handle Diameter: .25" (6,25 mm)

1025 [Sterilization Case]



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe

Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362

INNOMED