





Duncan Metatarsal Clamp

Designed by Gregory S. Duncan, DPM

Designed to be used on bones of the foot to stabilize an osteotomy or fracture in the corrected position for fixation through the opening in the top of the clamp

May also be used for open reduction internal fixation for hand or fibula procedures.



PRODUCT NO'S:	
1638 [Large]	Overall Length: 7" (17,8 cm) Clamp Pads: 1.3" x .625" (3,3 cm x 1,6 cm)
1638-25 [Medium]	Overall Length: 6.5" (16,5 cm) Clamp Pads: 1" x .5" (2,5 cm x 1,3 cm)
1638-50 [Small]	Overall Length: 6.25" (15,9 cm) Clamp Pads: .625" x .325" (1,6 cm x .8 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED 



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362