


Evans Modified Fukuda-type Retractors

Designed by Peter J. Evans, MD

*Designed to retract the humeral shaft posteriorly,
helping to expose the glenoid surface*

Center groove allows the shaft of a reamer to fit more posteriorly.



PRODUCT NO'S:

5180-N [Narrow]

Overall Length: 8.625"

Blade Width: 1" (25.4mm)

Blade Depth: 3.75"

5180-W [Wide]

Overall Length: 8.625"

Blade Width: 1.25" (31.7mm)

Blade Depth: 3.5"

MADE IN USA

ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362