

Evans Reverse Hohmann Retractor

Designed by Peter J. Evans, MD

Smaller size useful for retracting the deltoid superiorly or laterally, and also protecting the axillary nerve inferiorly while simultaneously exposing the glenoid



PRODUCT NO:
4547
Overall Length: 8.5" (21,6 cm)
Blade Width: Tapers from 30 mm to 18 mm
Blade Depth: 3" (7,6 cm)
Prong Width: 6 mm



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED 



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362