

# Gelbke Cobb Elevator with Suction

Designed by Martin K. Gelbke, MD

*Designed to be used during exposure of the posterior spine,  
as well as for pelvic and acetabular trauma cases*

Ultra hard titanium nitride coating helps to extend life by increasing surface hardness, prolonging sharpness, and resisting chemicals and corrosion.

**PRODUCT NO:**

3433

Overall Length: 12.75" (32,4 cm)

Cobb End Width: 18 mm (.7")

Shaft plus Head Length: 5.5" (15 cm)



ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

# INNOMED



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**1.800.548.2362**