



# George Semi-Circumferential Glenoid Retractor

Designed by Michael S. George, MD

*Designed to depress the humeral head and retract tissue away from the posterior half of the glenoid, helping to improve exposure for the preparation and placement of the glenoid component in total shoulder arthroplasty*

**PRODUCT NO:**

2435

Overall Length: 8" (20,3cm)

Blade Width: 2.125" (5,4cm)

**PROUDLY  
MADE  
IN THE  
USA**



ISO 9001:2008 • ISO 13485:2003

**FREE TRIAL ON MOST INSTRUMENTS**

Scan to  
Launch Our  
Website



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**INNOMED**

**1.800.548.2362**