



Glenosphere Component Retractor

Designed by Tim Seachris

Designed for use in total and reverse shoulder arthroplasty

Coated version helps to protect component surfaces.

PRODUCT NO:

5841 [Coated End]
Overall Length: 10.25" (26,7 cm)
Blade Width: .9375" (2,4 cm)

5841-01 [Uncoated End]
Overall Length: 10.25" (26,7 cm)
Blade Width: .9375" (2,4 cm)



INNOMED 5841-01 CE
(01)00840277125408 (10)12191

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362