


# Gunther Glenoid Retractor

Designed by Stephen B. Gunther, MD

*Ergonomic design helps to retract the humeral head posteriorly during glenoid exposure while avoiding reamer contact during shoulder replacement surgery*



**PRODUCT NO:**

1999

Overall Length: 11" (27,9 cm)

Neck Width: .625" (15,9 mm)

Prong Outside Width: 1" (25,4 mm)

Prong Inside Width: .625" (15,9 mm)



USA MADE

ISO 9001:2008 • ISO 13485:2003

**FREE TRIAL ON MOST INSTRUMENTS**

Scan to  
Launch Our  
Website



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**INNOMED**

**1.800.548.2362**