


Hawkins Shoulder Instruments

Designed by Richard J. Hawkins, MD

Designed to enhance exposure during shoulder arthroplasty procedures



| PRODUCT NO'S: | |
|---------------|---|
| 5090 | [Small Spreader w/Articular Arms] Overall Length: 6.25" (15,9 cm) Arm Depth: 2.25" (5,7 cm) Prong Width: 21 mm Prong Length: 16 mm |
| 5091 | [Large Spreader w/Articular Arms] Overall Length: 10.5" (26,7 cm) Arm Depth: 2.375" (6 cm) Prong Width: 23 mm Prong Length: 23 mm |
| 5092 | [Anterior Capsular Retractor] Overall Length: 11.25" (28,6 cm) Handle Length: 5.25" (13,3 cm) Blade Depth: 3.25" (8,3 cm) Blade Width: 19 mm |
| 5093 | [Small Pectoralis Retractor] Overall Length: 10.25" (26 cm) Handle Length: 5.25" (13,3 cm) Blade Depth: 2.5" (6,4 cm) Blade Width: 25 mm |
| 5094 | [Extra Small Pectoralis Retractor] Overall Length: 11" (27,9 cm) Handle Length: 5.25" (13,3 cm) Blade Depth: 1.5" (3,8 cm) Blade Width: 25 mm |
| 5095 | [Cobb Elevator] Overall Length: 11" (27,9 cm) Handle Length: 5.5" (14 cm) Blade Width: 19 mm |
| 5096 | [Humeral Head Retractor] Overall Length: 9" (22,9 cm) Blade Depth: 2.75" (7 cm) Blade Width: 37 mm |
| 5097 | [Anterior Glenoid Retractor] Overall Length: 11" (27,9 cm) Blade Depth: 2.75" (7 cm) Blade Width @ Fat Pad: 34 mm Blade Width @ Neck: 18 mm |
| 5098 | [Deltoid Retractor] Overall Length: 9.5" (24,1 cm) Blade Depth: 3.75" (9,5 cm) Blade Width @ Fat Pad: 45 mm Blade Width @ Neck: 32 mm |
| 5099 | [Modified Darrach Retractor] Overall Length: 10.75" (27,3 cm) Blade Width: 19 mm |



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to Launch Our Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362