


# Modified Hohmann Retractors

Handle is contoured to allow better leverage and visualization

Useful for retracting tissues at the margins of the joint. Can be passed over the margins of the joint and held in place with weights or by hand.

The OrthoLucent™ version is made of a strong, lightweight carbon fiber PEEK composite material, which is completely radiolucent, helps to prevent from marring component surfaces, and can be steam sterilized.

PRODUCT NO'S:
4535 [Narrow] Overall Length: 10" (25,4 cm) Blade Width: 14 mm
4535-R* [OrthoLucent™ Narrow] Overall Length: 10" (25,4 cm) Blade Width: 18 mm
4535-01 [Extra Deep Narrow] Overall Length: 11.625" (29,5 cm) Blade Width: 16,4 mm
4545 [Short-tipped Narrow] Designed by Carl DiRaimondo, MD Overall Length: 9.5" (24,1 cm) Blade Width: 14 mm
6595 [Wide] Overall Length: 10" (25,4 cm) Blade Width: 42,5 mm
6595-01 [Extra Deep Wide] Overall Length: 11.5" (29,2 cm) Blade Width: 42,5 mm



[www.innomed.net](http://www.innomed.net)

FREE TRIAL ON MOST INSTRUMENTS



ISO 13485:2016



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.