


# Horseshoe Shoulder Frame and Blade Assembly

*Designed to enhance exposure during shoulder arthroplasty procedures – once in place the frame and blades provide self-retaining exposure, helping to free up surgical assistants*



Horseshoe Frame



Horseshoe Blade - Wide




Horseshoe Blade - Narrow



Horseshoe Blade - Narrow, Bent



Horseshoe Blade - Small Narrow



Horseshoe Blade - Inferior



Delto Pectoral Approach



Supero Lateral Approach

PRODUCT NO'S:
2030-00 [Set] Set includes (1) Frame, (1) of Each Blade Style
<b>Also available individually:</b>
2030-01 [Horseshoe Frame] Overall Dimensions: 7" x 5" (17,8 cm x 12,7 cm) Frame Width: .7" (15 mm)
2030-02 [Blade - Wide] Blade Width: 22 mm Overall Length: 11" (27,9 cm)
2030-03 [Blade - Narrow] Blade Width: 14 mm Overall Length: 11" (27,9 cm)
2030-04 [Blade - Narrow, Bent] Blade Width: 14 mm Overall Length: 10" (25,4 cm) Handle Length: 4.5" (11,4 cm)
2030-05 [Blade - Small Narrow] Blade Width: 16 mm Blade Depth: 2" Overall Length: 8.5" (21,6 cm)
2030-06 [Blade - Inferior] Blade Width: Outside 34 mm, Inside 24 mm Overall Length: 11.5" (29,2 cm) Prong Length: 28 mm



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to Launch Our Website



**INNOMED** 



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

1.800.548.2362