

J.B. Redler Retractor

Designed by M.R. Redler, MD

*Uniquely balanced retractor for bone exposure
for a multitude of upper extremity procedures*

Double-angle design allows for ideal exposure with minimal effort to hold the retractor, while the assistant's hands are well out of the way of the exposure. The aperture in the base of the handle allows the retractor to be attached via a Penrose drain to the table for hands-free approach.

PRODUCT NO:

1645

Overall Length: 5" (12,7 cm)

MADE EXCLUSIVELY
FOR INNOMED IN
GERMANY



www.innomed.net

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

INNOMED 