

# Jana Anterior Hip Starter Rasp

Designed by Ajoy Jana, MD

*Designed to help begin preparation of the femoral canal prior to stem broaching – features a large offset handle and a striking plate*




- ▶ Strike plate if needed
- ▶ Blunt tip and gentle curve to avoid cortical penetration
- ▶ Long handle to give added leverage entering dense bone and easier use in large patients
- ▶ Optimized for use in direct anterior total hip, but can be used with any surgical approach

#### PRODUCT NO:

**3010**

Overall Length: 14.5" (36,8 cm)  
Rasp Shaft Total Length: 13" (33 cm)  
Length of Rasp Section: 4" (10,2 cm)  
Shaft Diameter: .375" (0,95 cm)  
Handle Length: 5.5" (24 cm)  
Striking Platform: 1.5" x 1.5" (3,8 x 3,8 cm)



[www.innomed.net](http://www.innomed.net)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

**INNOMED** 



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362