

# Joint, Calcaneal, Small Bone Compressor/Distractor

*Selection lever switches the mechanism from compression to distraction*

Simply squeeze the handle one time after direction selection to engage the mechanism. Two hole sizes for pin size selection. Calibrated ratchet handle.



**PRODUCT NO'S:**

4865-LS [Standard]  
Overall Length: 8.5" (21,6 cm)  
Holes For: .062" & .094" (1.6 & 2.5 mm) K-wire Pins



[www.innomed.net](http://www.innomed.net)

FREE TRIAL ON MOST INSTRUMENTS

**INNOMED** 

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362