


# Joint, Calcaneal and Small Bone Compressors

*Designed for compression in fracture and osteotomy procedures*

*Two hole sizes for ease of pin size selection:  
.062" (1.6 mm) & .094" (2.4 mm)*



PRODUCT NO'S:
4210-SC [Small] Overall Length: 6" (15,2 cm)
4210-XSC [Extra Small] Overall Length: 4.25" (10,8 cm)



[www.innomed.net](http://www.innomed.net)

FREE TRIAL ON MOST INSTRUMENTS

# INNOMED

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362