



Kaminsky OrthoLucent™ Browne-type Deltoid Retractors

Designed by Sean B. Kaminsky, MD

Used for the Delto-Pectoral approach, and can also remain in place for fracture reduction, plate positioning, and screw/wire/drill location confirmation

Used for acromioplasty, rotator cuff repair, and fracture fixation. Contours the humeral head with deltoid retraction allowing extensive exposure. Helps to reduce operative time, assist in fracture reduction, and maintain hardware position without the frequent need for retractor removal and reintroduction.

The OrthoLucent™ carbon fiber PEI composite material is strong, lightweight, completely radiolucent, helps to prevent from marring component surfaces, and can be steam sterilized



PRODUCT NO'S:
1670-01R [Small] Blade Width: 4,5 cm Overall Length: 10.5" (26,7 cm)
1670-02R [Large] Blade Width: 5,4 cm Overall Length: 10.5" (26,7 cm)

MADE EXCLUSIVELY
FOR INNOMED IN
SWITZERLAND

ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Scan to
Launch Our
Website



INNOMED

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362

© 2015 Innomed, Inc.