



Kopplin Osteophyte Rongeur for Direct Anterior THA

Designed by Matthew Kopplin, MD


Designed to help remove osteophytes around the acetabulum in anterior THA

Helps to allow clean sharp cuts with better control. Thin flat tip helps to pass along the bone easier.



PRODUCT NO:
1771-04
Overall Length: 13" (35 cm)
Jaw Width: .375" (0.95 cm)
Jaw Length Bottom: .75" (1.9 cm)
Jaw Length Top: .625" (1.6 cm)

MADE EXCLUSIVELY
FOR INNOMED IN
GERMANY



www.innomed.net



FREE TRIAL ON MOST INSTRUMENTS

INNOMED

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2026 Innomed, Inc.