

# Large T-Handle Fixed Drivers

Large easy grip soft silicone handled drivers help provide a sturdy non-slip grip




Zimmer Hall Quick-connect

Zimmer Hall Reverse Quick-connect

Hudson Quick-connect

The two standard Quick-connect models release by pulling the collar backward, while the Reverse Quick-connect model is designed to have the collar be pushed forward for release.



PRODUCT NO'S:	
8248	[Zimmer Hall Quick-connect] Overall Length: 5.75" (15,6 cm) Handle Width: 4.625" (11,6 cm)
8248-01	[Reverse Quick-connect Zimmer Hall] Overall Length: 5.75" (15,6 cm) Handle Width: 4.625" (11,6 cm)
8249	[Hudson Quick-connect] Overall Length: 6.75" (17,1 cm) Overall Length with Pin In Handle: 11.5" (29,2 cm)

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to Launch Our Website



**INNOMED**



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

1.800.548.2362