

# Lighted Yankaur Suction Device

Designed by Adolph V. Lombardi Jr., MD

*Designed to help provide effective suction with the addition of a light source for enhanced visualization*

**NOTE: The Disposable LED Light Source is currently unavailable. The device can still be used with a fiber optic light cable.**




- ▶ Device comes with one (1) Disposable LED Light Source (#8010-01)
- ▶ Can also be attached to a fiber optic light cable with ACMI (female) connector
- ▶ The handle is made of Delrin
- ▶ Entire device is steam sterilizable

## Disposable LED Light Source



<b>PRODUCT NO'S:</b>
<b>PACKAGE OF 1:</b>
8010-01 [Disposable LED Light Source] Overall Length: 2.5" (6,4 cm) Diameter: 1" (2,54 cm)
<b>PACKAGE OF 10:</b>
8010-10 [Disposable LED Light Source]

<b>PRODUCT NO'S:</b>
8016-L-01 [Kit – Device with Disposable LED Light Source] Overall Length: 11.75" (29,8 cm) Handle Length: 3.93" (10 cm) Handle Width: .86" (2,2 cm) Suction Tube Diameter: .25" (6,35 mm)
<b>Also available individually:</b>
8016-L [Suction Device Only] Overall Length: 11.75" (29,8 cm) Handle Length: 3.93" (10 cm) Handle Width: .86" (2,2 cm) Suction Tube Diameter: .25" (6,35 mm)



ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

# INNOMED



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**1.800.548.2362**

© 2021 Innomed, Inc.