

Long Handled Fukuda Type Retractors

Design modification of an original design by Evan Flatow, MD & Louis Bigliani, MD

Designed to retract the humeral shaft posteriorly, helping to expose the entire glenoid surface, these retractors feature a longer handle for use when added leverage and depth is needed for use in large patients



PRODUCT NO'S:	
1930-H	[Narrow with Long Handle] Blade Width: 32 mm Opening: 25 x 40 mm Overall Length: 9.35" (23,7 cm)
1940-H	[Wide with Long Handle] Blade Width: 38 mm Opening: 32 x 40 mm Overall Length: 9.35" (23,7 cm)



www.innomed.net



FREE TRIAL ON MOST INSTRUMENTS

INNOMED 

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2026 Innomed, Inc.