

# Lower Extremity Leg Positioner


Designed by Ronald Romanelli, MD

*Used to support knee and leg during surgery, and can be used for casting*

- ▶ Utilized for rodding of femurs or tibias
- ▶ Also useful for knee surgery and closures
- ▶ Very supportive, distributes stresses on leg, used instead of bolsters
- ▶ Supplied with one autoclavable silicone pad
- ▶ Aluminum positioner is radiolucent and gas or steam sterilizable



<b>PRODUCT NO'S:</b>
2745
Dimensions: 5.5" H x 9.5" L x 9.25" W (12,7 cm x 24,1 cm x 23,5 cm)
<b>Replacement Parts:</b>
2760-P [Silicone Pad]



[www.innomed.net](http://www.innomed.net)

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

**INNOMED** 



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362