

Meyerding Type Retractors with Contour Grip

Designed for soft tissue retraction, the contour grip allows for a better grip and less fatigue

Features a non-glare finish.



50 x 16 mm Blade
#6241-01



75 x 15 mm Blade
#6242-01



75 x 25 mm Blade
#6243-01

PRODUCT NO'S:	
6241-01	[50 x 16 mm] Overall Length: 8.5" (21,6 cm) Blade Width: 16 mm Blade Depth: 50 mm
6242-01	[75 x 15 mm] Overall Length: 8.75" (22,2 cm) Blade Width: 15 mm Blade Depth: 75 mm
6243-01	[75 x 25 mm] Overall Length: 8.75" (22,2 cm) Blade Width: 25 mm Blade Depth: 75 mm

MADE EXCLUSIVELY
FOR INNOMED IN
GERMANY



– See other side for alternate handle style –

www.innomed.net



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED 



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2026 Innomed, Inc.

Meyerding Type Retractors with Ergonomic Handle

Designed for soft tissue retraction, the ergonomic handle allows for a better grip and less fatigue

Non-glare finish featured on the metal retractor parts.



PRODUCT NO:

6242 [75 x 15 mm]
Overall Length: 9" (22.9 cm)
Blade Width: 15 mm
Blade Depth: 75 mm

6243 [75 x 25 mm]
Overall Length: 9" (22.9 cm)
Blade Width: 25 mm
Blade Depth: 75 mm

– See other side for alternate handle style –

www.innomed.net



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2026 Innomed, Inc.