


Micro Curettes

Four cup sizes, straight or 45° angled-end shaft



Ultra hard titanium nitride coating helps to extend curette life by increasing surface hardness, prolonging sharpness, and resisting chemicals and corrosion.

PRODUCT NO'S:

Straight Micro Curettes	Angled Micro Curettes
Overall Length: 9.75" Shaft Length: 4.5"	Overall Length: 9.75" Shaft Length: 4.5"
4242 <i>Cup Size 2</i>	4242-01 <i>Cup Size 2</i>
4240 <i>Cup Size 1</i>	4240-01 <i>Cup Size 1</i>
4244 <i>Cup Size 4/0</i>	4244-01 <i>Cup Size 4/0</i>
4246 <i>Cup Size 6/0</i>	4246-01 <i>Cup Size 6/0</i>

**PROUDLY
MADE
IN THE
USA**

ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362