

Modified Fukuda-type Retractor with Reamer Slot

Designed by Richard J. Miller, MD

Center cutout slot allows the shaft of a reamer to fit more posteriorly

Used to retract the humeral shaft posteriorly and helping to expose the entire glenoid surface.
Available in small and large sizes.



PRODUCT NO'S:

1898 [Small]

Overall Length: 7.25"

Blade Width: 32mm

Opening: 25x40mm

1899 [Large]

Overall Length: 7.25"

Blade Width: 38mm

Opening: 32x40mm

**PROUDLY
MADE
IN THE
USA**

ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362