


Modular Impactor Set

Makes multiple impactor heads easily visible and available


Designed to have available to the operating surgeon multiple types of impactors utilizing one handle. The rack uses less space and allows the surgeon to quickly see the designs available. The impactors are supplied with stainless steel tips for bone and delrin tips which can be used against an implant for slight placement adjustments.



PRODUCT NO:	
5370	[Complete Set]
Included In Set / Also Available Individually:	
5370-01	[Rectangular Impactor Tip 11mm x 4 mm Steel]
5370-02	[Oval Impactor Tip 13 mm x 8 mm Steel]
5370-03	[Crescent Impactor Tip 12 mm x 5 mm Steel]
5370-04	[Square Impactor Tip 9 mm x 9 mm Steel]
5370-05	[Round Impactor Tip 15 mm Steel]
5370-06	[Round Impactor Tip 12 mm Steel]
5370-07	[Round Impactor Tip 9 mm Steel]
5370-19	[Impactor Set Base] Base Diameter: 3.5" (8,9 cm)
5370-D1	[Rectangular Impactor Tip 11mm 4 mm Delrin]
5370-D2	[Oval Impactor Tip 13 mm x 8 mm Delrin]
5370-D3	[Crescent Impactor Tip 12 mm x 5 mm Delrin]
5370-H	[Modular Impactor Handle] Overall Length: 8" (20,3 cm) Grip Length: 4.5" (11,4 cm)



Stainless Steel
Impactor Sizes



Delrin
Impactor Sizes



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362

© 2021 Innomed, Inc.