


# Mueller-type Femoral Neck Elevator

Extra Deep modified by Tom Eickmann, MD

*Designed to elevate the proximal femur*



Standard  
#3415

Rounded Tips  
#3415-02



Rounded tips help prevent the tips from digging into bone

Extra Deep  
#3418

| PRODUCT NO:                   |   |
|-------------------------------|---|
| <b>3415</b> [Standard]        | Overall Length: 13.5" (34,3 cm)<br>Handle Length: 6.5" (16,5 cm)<br>Blade Width at Widest: 25 mm  |
| <b>3415-02</b> [Rounded Tips] | Overall Length: 13.375" (34 cm)<br>Handle Length: 6.5" (16,5 cm)<br>Blade Width at Widest: 25 mm  |
| <b>3418</b> [Extra Deep]      | Overall Length: 15.25" (38,8 cm)<br>Handle Length: 6.5" (16,5 cm)<br>Blade Width at Widest: 25 mm |



[www.innomed.net](http://www.innomed.net)



# INNOMED

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362