



# Offset Handle Hip Retractor

Designed by Gina Bart, PAC

*The 90° offset handle helps move the assistant's hand out of the way, allowing better access to the surgical site*



During posterior THR, the instrument allows for one assistant to hold the leg while retracting the abductors, thereby maximizing both a clearer view and better physical access to the surgical site.

**PRODUCT NO:**

6207-L [Left]	6207-R [Right]
Overall Length: 13"	Overall Length: 13"
Handle Offset: 8"	Handle Offset: 8"
Blade Width: 15mm	Blade Width: 15mm



**FREE TRIAL ON MOST INSTRUMENTS**



103 Estus Drive, Savannah, GA 31404  
[www.innomed.net](http://www.innomed.net) [info@innomed.net](mailto:info@innomed.net)

912.236.0000 Phone  
 912.236.7766 Fax

© 2012 Innomed, Inc.

**Innomed-Europe** Tel. +41 41 740 67 73  
 Fax +41 41 740 67 74



**1.800.548.2362**