

# Offset Punches

*Helps in the removal of hip stems*

Used to help remove a hip prosthesis stem via a window in the shaft of the femur. Two sizes of offsets allow the punches to be used to tap on a distal portion of the hip stem, after a window has been made in the femur below the tip of the stem.

#### PRODUCT NO'S:

5125-01 [Small Offset]  
Overall Length: 11" (27,9 cm)  
Punch End Offset: .5" (1,3 cm)  
Punch End Diameter: 7 mm

5125-02 [Large Offset]  
Overall Length: 11" (27,9 cm)  
Punch End Offset: 1.25" (3,2 cm)  
Punch End Diameter: 7 mm



[www.innomed.net](http://www.innomed.net)

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

# INNOMED



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362