


# OrthoLucent™ Hohmann Retractors

Carbon fiber PEEK composite material is lightweight, completely radiolucent and can be steam sterilized. Also helps to prevent from marring component surfaces.



Modified Narrow Hohmann  
#4535-R

Modified Blunt Hohmann  
#4550-R

Standard Hohmann  
#4558-R

Narrow Bent Hohmann  
#7110-R

PRODUCT NO'S:
4535-R [Modified Narrow] Blade Width: 18 mm Overall Length: 10" (25,4 cm)
4550-R [Modified Blunt] Blade Width at Widest: 24,5 mm Overall Length: 10.75" (27,3 cm)
4558-R [Standard] Blade Width: 16 mm Overall Length: 9.46" (24 cm)
7110-R [Narrow Bent] Overall Length: 9.75" (23,8 cm) Handle Length: 7" (17,8 cm) Blade Width: 19 mm Depth from Bend: 4.75" (12,1 cm)

MADE EXCLUSIVELY  
FOR INNOMED IN  
SWITZERLAND

[www.innomed.net](http://www.innomed.net)

FREE TRIAL ON MOST INSTRUMENTS

# INNOMED

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362