


Ortho Self-Retaining Retractor with Pin Guides

Designed by Sean Dunn, DPM

Designed for small joint use with pin guides that are set back to allow either direct distraction or distraction with pins

Uses include:

- ▶ Osteotomy distraction (such as the Evans or Cotton in the foot)
- ▶ Joint distraction for arthodesis or lengthening applications
- ▶ Fracture distraction



- ▶ Parallel pin guides allow pins up to 2 mm
- ▶ Serrated outside blades extend .4" (1 cm) beyond end of pin guides

PRODUCT NO:

1842-02

Overall Length: 6.5" (16,5 cm)

Blade Width: 7 mm

Blade Extension (beyond guides): .4" (1 cm)

Blade Thickness: 1.68 mm

Pin Guide Length: 1.25" (3,2 cm)

Pin Guide Internal Diameter: .085" (2,1 mm)



Scan to
Launch Our
Website



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362