


# Proximal Tibia Resection Shark Tooth Clamp


Designed by Shara Diers, PA-C

*Designed to help grasp and remove the cut proximal portion of the tibia during total and uni knee arthroplasty, the small version is also useful in total ankle surgery*

The top looping jaw of the standard size clamp features a sharp forked end for grasping, while the small version features a single pointed end.



Standard  
#3651




Can also be useful in  
total ankle surgery

Small  
#3651-S



PRODUCT NO'S:
<b>3651 [Standard]</b> Overall Length: 7" (17,8 cm) Platform Width: .79" (2 cm) Bottom Jaw Length: 1.5" (3,8 cm)
<b>3651-S [Small]</b> Overall Length: 7" (17,7 cm) Platform Width: .25" (0,63 cm) Bottom Jaw Length: .75" (1,9 cm)

MADE EXCLUSIVELY  
FOR INNOMED IN  
GERMANY



[www.innomed.net](http://www.innomed.net)



FREE TRIAL ON MOST INSTRUMENTS

# INNOMED

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362

© 2026 Innomed, Inc.