


# Radiolucent Self-Retaining Hip Surgery Retractor Frame & Blade Set

*A radiolucent locking frame and blade set for use during total hip arthroplasty*

Compatible with all retractors of the Innomed Self-Retaining Hip Surgery Retractor System.



Radiolucent  
Locking Frame  
#7325-01

*Completely radiolucent with anodized  
aluminum handles and delrin blades*

2" Radiolucent  
Blade  
#7449-02R

3" Radiolucent  
Blade  
#7449-03R

4" Radiolucent  
Blade  
#7449-04R

<b>PRODUCT NO'S:</b>
<b>7325-00</b> [Frame & 3 Blade Sizes]
<b>Also available Individually:</b>
<b>7325-01</b> [Radiolucent Locking Frame] 9.75" x 7" (24,8 cm x 217,8 cm)
<b>7449-02R</b> [2" Radiolucent Blade] Blade Depth: 2" (5,1 cm) Handle Length: 6" (15,2 cm)
<b>7449-03R</b> [3" Radiolucent Blade] Blade Depth: 3" (7,6 cm) Handle Length: 6" (15,2 cm)
<b>7449-04R</b> [2" Radiolucent Blade] Blade Depth: 4" (10,2 cm) Handle Length: 6" (15,2 cm)



[www.innomed.net](http://www.innomed.net)



# INNOMED

ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
info@innomed-europe.com

**1.800.548.2362**