

Robb Cement Curette

Designed by William Robb, MD

Designed to help remove cement around a knee or hip prosthesis

PRODUCT NO:

5635

Overall Length: 8" (20,3 cm)

Freer End: 5 mm

Cup End: 10 mm

Made of Delrin



www.innomed.net

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.