

# Rogozinski Reverse Angle Retractors

Designed by Chaim Rogozinski, MD

*Designed to be self-leveling, helping to maintain the body of the retractor on the patient for soft tissue retraction and out of the surgeons field, with finger loops designed for use with either hand*

Designed for spine but can be used for other surgeries as well.

## PRODUCT NO'S:


4272 [Large]  
 Overall Length: 9" (22,9 cm)  
 Length to Bend: 8.5" (21,6 cm)  
 Depth: 4.25" (10,8 cm)

4273 [Medium]  
 Overall Length: 8" (20,3 cm)  
 Length to Bend: 8" (20,3 cm)  
 Depth: 3" (7,6 cm)

4274 [Small]  
 Overall Length: 8" (20,3 cm)  
 Length to Bend: 8" (20,3 cm)  
 Depth: 1.75" (4,4 cm)



*Typical position of most retractor handles*



*Design allows handle to fall back towards the patient and out of the surgeon's way*

ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

Scan to  
 Launch Our  
 Website



**INNOMED**



103 Estus Drive, Savannah, GA 31404  
 www.innomed.net info@innomed.net

912.236.0000 Phone  
 912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
 Fax +41 41 740 67 71

**1.800.548.2362**

© 2018 Innomed, Inc.