

Rogozinski Capsule Dilator

Designed by Chaim Rogozinski, MD






A reverse plier mechanism with a pointed tip used to puncture and help pry open the capsule and to dilate an arthroscopic portal

PRODUCT NO:

5186

Overall Length: 6.875" (17,5 cm)
Jaw Length: 3.5" (8,9 cm)

MADE EXCLUSIVELY
FOR INNOMED IN
GERMANY



www.innomed.net



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362